

Capp *Harmony*

SEROLOGICAL PIPETTES

Capp *Harmony* pipettes are made from transparent medical-grade polystyrene sterilized by gamma irradiation. Each lot is tested for endotoxins (C.5 Eu/MI) (ISO TO993-11) and Cytotoxicity (ISO 10993-5) and are guaranteed non-haemolytic (ISO 10993-4). Ideal for use within tissue culture applications.

Features:

- Designed for accurate and rapid fluid dispensing, these polystyrene pipettes are sterile, pyrogen-free and cotton-plugged for a variety of liquid handling applications
- The sharp, black, highly legible permanent graduations prevent ambiguity and guarantee precise reading and accurate dispensing to +1-2%.
- Pipettes are supplied individually wrapped in fibre-free colour-coded paper and plastic combination that is easy to open
- Each pipette comes with an effective filter barrier to prevent overflow of liquid into pipette
- Reverse graduations show the volume removed or remaining on all pipettes larger than 1mL.
- Generous negative graduations provide additional pipetting volume.
- The uniform mouth piece fits all major pipettors and is colour-coded for ease of identification



Ordering Information

Cat. No.	Description	Volume	Package	Colour
SP-1-C	Serological pipette	1mL	500 pcs.	Yellow
SP-2-C	Serological pipette	2mL	500 pcs.	Green
SP-5-C	Serological pipette	5mL	200 pcs.	Cyan
SP-10-C	Serological pipette	10mL	200 pcs.	Orange
SP-25-C	Serological pipette	25mL	150 pcs.	Blue
SP-50-C	Serological pipette	50mL	100 pcs	Purple
SP-100-C	Serological pipette	100mL	50 pcs.	Pink

Capp *Harmony*

RETIPS

Capp offers a wide range of **Capp *Harmony*** retips, featuring universal compatibility with the majority of repeater pipettes available on the market.

The concept of the retips is based on the positive displacement principle, where the piston of the retip is in direct contact with the liquid being processed, minimizing the dead volume leading to increased inaccuracy and imprecision margins.

Capp *Harmony* retips are therefore ideal to assure very reliable results on your liquid handling protocol. Moreover, the disposable pistons of **Capp *Harmony*** retips are perfectly suitable for applications involving viscous liquids and for protocols requiring an increased prevention of cross contamination, such as pipetting of infectious and volatile liquids.

Capp *Harmony* retips along with Capp repeater pipettes - **Capp *Rhythm*** or **Capp *R10*** - are the optimal tool for any multidispensing application. Available in 11 sizes (8 for **Capp *Rhythm***, 11 for **Capp *R10*** repeater) they come bulk or in an individual sterile packaging.



Ordering Information

Capp *Harmony* Syringes

Cat. No. Description

RS-01-C	Syringes 0,05 mL, 100 pcs.
RS-11-C	Syringes Pro 0,1 mL, 100 pcs.
RS-12-C	Syringes Pro 0,2 mL, 100 pcs.
RS-02-C	Syringes Pro 0,5 mL, 100 pcs.
RS-13-C	Syringes Pro 1,0 mL, 100 pcs.
RS-03-C	Syringes 1,25 mL, 100 pcs.
RS-04-C	Syringes Pro 2,5 mL, 100 pcs.
RS-05-C	Syringes Pro 5,0 mL, 100 pcs.
RS-14-C	Syringes Pro 10,0 mL, 100 pcs.
RS-06-C	Syringes 12,5 mL, 100 pcs.
RS-07-C	Syringes Pro 25,0 mL, 25 pcs.
RS-08-C	Syringes Pro 50,0 mL, 25 pcs.

Capp *Harmony* Syringes, Sterile, Individually Packed

Cat. No. Description

RS-01S-C	Syringes 0,05 mL, sterile, individually packed, 100 pcs.
RS-11S-C	Syringes Pro 0,1 mL, sterile, individually packed, 100 pcs.
RS-12S-C	Syringes Pro 0,2 mL, sterile, individually packed, 100 pcs.
RS-02S-C	Syringes Pro 0,5 mL, sterile, individually packed, 100 pcs.
RS-13S-C	Syringes Pro 1,0 mL, sterile, individually packed, 100 pcs.
RS-03S-C	Syringes 1,25 mL, sterile, individually packed, 100 pcs.
RS-04S-C	Syringes Pro 2,5 mL, sterile, individually packed, 100 pcs.
RS-05S-C	Syringes Pro 5,0 mL, sterile, individually packed, 100 pcs.
RS-14S-C	Syringes Pro 10,0 mL, sterile, individually packed, 100 pcs.
RS-06S-C	Syringes 12,5 mL, sterile, individually packed, 100 pcs.
RS-07S-C	Syringes Pro 25,0 mL, sterile, individually packed, 25 pcs
RS-08S-C	Syringes Pro 50,0 mL, sterile, individually packed, 25 pcs.

Capp Harmony

COMPATIBILITY CHART

For syringes' compatibility with other pipette brands, please contact info@capp.dk



RS-01-C

RS-01S-C

Capp Harmony syringes 0.05mL
100 pcs



RS-11-C

RS-11S-C

Capp Harmony syringes 0.1mL
100 pcs



RS-12-C

RS-12S-C

Capp Harmony syringes 0.2mL
100 pcs



RS-02-C

RS-02S-C

Capp Harmony syringes 0.5mL
100 pcs



RS-13-C

RS-13S-C

Capp Harmony syringes 1.0mL
100 pcs



RS-03-C

RS-03S-C

Capp Harmony syringes 1.25mL
100 pcs

Capp Rhythm Repeater	Fully compatible	Not compatible	Not compatible	Fully compatible	Not compatible	Fully compatible
Capp R10 Repeater	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Biohit Mechanical Stepper	Fully compatible	Not compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible
eLINE Lite Dispenser	Not compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Not compatible
eLINE Pro Dispenser	Not compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Not compatible
Fisherbrand	Fully compatible	Not compatible	Not compatible	Fully compatible	Not compatible	Fully compatible
HTL MINILAB 201	Fully compatible	Not compatible	Not compatible	Fully compatible	Not compatible	Fully compatible
Brand HandyStep	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Brand HandyStep eletronic	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Not compatible	Fully compatible
Gilson Distriman / Repetman	Not compatible	Not compatible	Not compatible	Not compatible	Not compatible	Fully compatible
Raining Auto Rep S	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Raining Auto Rep E	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Not compatible	Fully compatible

Capp Harmony

COMPATIBILITY CHART

For syringes' compatibility with other pipette brands, please contact info@capp.dk



RS-04-C RS-04S-C Capp Harmony syringes 2.5mL 100 pcs	RS-05-C RS-05S-C Capp Harmony syringes 5mL 100 pcs	RS-14-C RS-14S-C Capp Harmony syringes 10mL 100 pcs	RS-06-C RS-06S-C Capp Harmony syringes 12.5mL 100 pcs	RS-07-C RS-07S-C Capp Harmony syringes 25mL 25 pcs	RS-08-C RS-08S-C Capp Harmony syringes 50mL 25 pcs
Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible
Not compatible	Not compatible	Not compatible	Fully compatible	Not compatible	Not compatible
Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible