

Part No.	Kinematic Viscosity mm ² /s (cSt)	
	20.00°C	25.00°C
	68.00°F	77.00°F
ALK-VIS-RT5-600	5.743	5.278
ALK-VIS-RT10-600	11.81	10.79
ALK-VIS-RT50-600	58.53	53.42
ALK-VIS-RT75-600	88.21	80.38
ALK-VIS-RT100-600	115.2	104.3
ALK-VIS-RT250-600	289.8	263.9
ALK-VIS-RT350-600	410.4	373.4
ALK-VIS-RT500-600	579.6	527.0
ALK-VIS-RT1K-600	1157	1051
ALK-VIS-RT5K-600	5840	5296
ALK-VIS-RT10K-600	11341	10273
ALK-VIS-RT12K-600	13745	12446
ALK-VIS-RT30K-600	33389	30244
ALK-VIS-RT60K-600	64845	58742
ALK-VIS-RT100K-600	108039	97858

Part No.	Dynamic Viscosity mPa's (cP)	
	20.00°C	25.00°C
	68.00°F	77.00°F
ALK-VIS-RT5-600	5.278	4.826
ALK-VIS-RT10-600	11.09	10.08
ALK-VIS-RT50-600	55.14	50.08
ALK-VIS-RT75-600	83.45	75.67
ALK-VIS-RT100-600	111.5	100.5
ALK-VIS-RT250-600	275.2	249.4
ALK-VIS-RT350-600	390.1	353.2
ALK-VIS-RT500-600	552.6	500.0
ALK-VIS-RT1K-600	1107	1001
ALK-VIS-RT5K-600	5651	5101
ALK-VIS-RT10K-600	11034	9949
ALK-VIS-RT12K-600	13399	12076
ALK-VIS-RT30K-600	32541	29343
ALK-VIS-RT60K-600	63217	57003
ALK-VIS-RT100K-600	105327	94961

Part No.	Density g/mL	
	20.00°C	25.00°C
	68.00°F	77.00°F
ALK-VIS-RT5-600	0.9190	0.9143
ALK-VIS-RT10-600	0.9389	0.9343
ALK-VIS-RT50-600	0.9421	0.9375
ALK-VIS-RT75-600	0.9460	0.9414
ALK-VIS-RT100-600	0.9678	0.9633
ALK-VIS-RT250-600	0.9495	0.9449
ALK-VIS-RT350-600	0.9505	0.9459
ALK-VIS-RT500-600	0.9534	0.9488
ALK-VIS-RT1K-600	0.9567	0.9521
ALK-VIS-RT5K-600	0.9677	0.9632
ALK-VIS-RT10K-600	0.9729	0.9685
ALK-VIS-RT12K-600	0.9748	0.9703
ALK-VIS-RT30K-600	0.9746	0.9702
ALK-VIS-RT60K-600	0.9749	0.9704
ALK-VIS-RT100K-600	0.9749	0.9704